HORIBA

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HORIBA

65□□-10C COMBINATION TYPE ION SELECTIVE ELECTRODE INSTRUCTION MANUAL

The combination type ion selective electrode is a 3-in-1 electrode that combines a tip type ion electrode and reference electrode. Because you do not need to deal with two electrodes (the ion electrode and the reference electrode), this is an extremely handy way to take ion concentration measurements. The ion electrode section uses a tip type electrode in a convenient structure that allows easy removal and replacement. Using combinations of the various types of tips, you can make many types of ion measurements with this single electrode. The reference electrode is a double junction structure with free selection of internal solution. This has the merit of minimizing contamination of the sample liquid from internal solution outflow, and minimizing electrode potential disturbances. The moveable sleeve type liquid junction ensures there is no blockage of the internal solution, and minimizes liquid junction potential.

Be sure to read these instructions carefully before using the equipment to ensure correct usage of the ion electrode.

1. Contents of the Package

- Electrode main unit 1 (A rubber cap is attached to the tip type ion electrode section at shipment.)
- package - Electrode tips 1 (The number of tips is listed in the individual specifications table below.)
- Syringe 1 (1 mL)
- Tip Type Electrode Instruction Manual 1 (Contains an explanation of the standard liquid adjustment method, detection characteristic diagrams, etc.)
- Instruction manual 1

2. Specifications

● Common specifications

- Temperature range : 0 to 50 degrees C.
- Internal electrode : Ag/AgCl (silver/silver chloride) electrode
- Internal solution (Inner tube) : 3.33 mol/L KCl gel solution (Cat. No. 330)
(Outer tube) : Differs depending on the type of ion being measured. (Listed in the individual specifications table).
- Liquid junction (Inner tube) : Multi-aperture high density polyethylene resin
(Outer tube) : Moveable sleeve
- Material of the parts that come in contact with the liquid : Poly vinyl chloride resin, polysulfone resin
*Cannot be used with samples that contain organic solvents (such as acetone, tetrahydrofuran, etc.)

● Individual specifications

Model Ion Type	Accessory Tip Type Electrode Model (Number of Tips)	① Concentration Range ② pH Range* ¹ ③ Slope (mV/activity decade)* ²	Interfering Ions Selectivity Constant* ³	Reference Electrode (Outer Tube) Type of Internal Solution* ⁴
#6560-10C Cl ⁻	#7660 (1 piece)	① 0.4~35,000mg/L Cl ⁻ (1~10 ⁻⁵ mol/L Cl ⁻) ② pH3~11 (350mg/L Cl ⁻) ③ -59mV	NO ₃ ⁻ , F ⁻ , HCO ₃ ⁻ , SO ₄ ²⁻ , PO ₄ ²⁻ =1000 以上 SCN ⁻ =0.3 MnO ₄ ⁻ =0.1 Br ⁻ =0.03 S ₂ O ₃ ²⁻ , S ²⁻ , I ⁻ , Ag ⁺ , Hg ²⁺ =Not possible (in 10 ⁻³ mol/L Cl ⁻)	1mol/L KNO ₃ (100 g/L KNO ₃)
#6561-10C F ⁻	#7661 (1 piece)	① 0.02~19,000mg/L F ⁻ (1~10 ⁻⁶ mol/L F ⁻) ② pH4~10 (20mg/L F ⁻) ③ -59mV	OH ⁻ =10 (within the concentration range)	3.33mol/L KCl (250 g/L KCl)
#6581-10C NO ₃ ⁻	#7681 (2 pieces)	① 0.62~62,000mg/L NO ₃ ⁻ (1~10 ⁻⁵ mol/L NO ₃ ⁻) ② pH3~7 (62mg/L NO ₃ ⁻) ③ -55mV	SO ₄ ²⁻ =1000 以上 CH ₃ COO ⁻ =300 F ⁻ =200 Cl ⁻ =40 NO ₂ ⁻ =3 I ⁻ =0.1 ClO ₄ ⁻ =0.03 Br ⁻ =2 (in 10 ⁻³ mol/L NO ₃ ⁻)	1mol/L KCl (75 g/L KNO ₃)
#6582-10C K ⁺	#7682 (2 pieces)	① 0.04~39,000mg/L K ⁺ (1~10 ⁻⁶ mol/L K ⁺) ② pH5~11 (3.9mg/L K ⁺) ③ +58mV	Li ⁺ , Na ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ =1000 以上 NH ₄ ⁺ =70 Cs ⁺ =3 Rb ⁺ =0.4 (in 10 ⁻⁴ mol/L K ⁺)	3.33mol/L NaCl (200 g/L KNO ₃)
#6583-10C Ca ²⁺	#7683 (2 pieces)	① 0.4~40,080mg/L Ca ²⁺ (1~10 ⁻⁵ mol/L Ca ²⁺) ② pH5~11 (4mg/L Ca ²⁺) ③ +29mV	Na ⁺ , K ⁺ , Ba ²⁺ , NH ₄ ⁺ =1000 以上 Mg ²⁺ =1000 Mn ²⁺ =500 Co ²⁺ =350 Ni ²⁺ , Cu ²⁺ =70 Sr ²⁺ =50 Fe ²⁺ , Zn ²⁺ =1 Fe ³⁺ =0.1 (in 10 ⁻⁴ mol/L Ca ²⁺)	3.33mol/L KCl (250 g/L KCl)

- (* 1) This is the allowed pH range for the measured ion concentration value given inside the brackets. (The pH range varies slightly depending on the ion concentration).
- (* 2) This is the change in the electrode potential when the ion concentration becomes 10 times denser (at 25 degrees C.)
- (* 3) When the concentration of coexisting ions exceeds the selectivity constant for the measured ion concentration, they will affect the electrode potential and make measurement impossible. (The effects of the other ions becomes smaller as K becomes larger, or when the ion concentration in the sample liquid increases.)
- (* 4) The type of internal solution that can be used is not limited to the listed solutions. Other solutions may be used if the outflow does not disturb the electrode potential, and does not generate any electric potential differences between liquids.

3. The Electrode


Figure 1. The combination type ion electrode

4. Preparations

- ① Manufacturing the internal solution
A commercially available guaranteed reagent grade reagent is used for manufacturing the internal solution. Measure the quantity of reagent given in the brackets in the internal solution column of the individual specifications table. Put it in a 1 liter flask and add pure (deionized) water to fill the flask. (An accuracy level of ± 2 g is adequate when measuring the reagent.)
- ② Manufacturing the standard solution
The procedure is described in the instruction manual for the tip type ion electrode.

5. Before Usage

- ① Remove the protective cap.
 - ② Remove the rubber stopper from the internal solution filling opening. With the electrode inverted, use the syringe to remove the internal solution.
 - ③ Use the syringe to fill the tube with the internal solution manufactured in step 4.(1).
 - ④ Remove the protective pipe. (Verify that the rubber cap is attached to the tip type ion electrode section.)
 - ⑤ If any crystals from the internal solution stick to the liquid junction, soak the crystals with pure (deionized) water until they dissolve, then wipe away any water droplets with gauze or tissue paper.
 - ⑥ Stand the rubber cap end of the electrode on a desk or other horizontal surface and press it down two or three times until the internal solution starts to come out. (Figure 5)
- Note:**
1. Check that the rubber stopper for the internal solution filling opening has been removed.
 2. Make sure that the rubber cap is attached for the tip type ion electrode section. If this operation is performed with the tip electrode attached, it may cause damage to the electrode.
- ⑦ Wash the electrode well with pure (deionized) water. Wipe away any water droplets with gauze or tissue paper.
 - ⑧ Remove the rubber stopper. Insert the tip type ion electrode and attach the protective pipe.
 - ⑨ Condition the electrode under the following conditions to prepare it for measurement. The solid state membrane electrodes #6560 (Cl^-) and #6561 (F^-) do not require any special conditioning. Before starting the measurements, soak the electrode in the high concentration standard liquid for approximately 3 minutes.

Type of Electrode	Conditioning Liquid	Soak Time
#6581 (NO_3^-)	1 mol/L KNO_3 solution (100g/L KNO_3)	Approx. 1Hr
#6582 (K^+)	0.1 mol/L KCl solution (75g/L KCl)	Approx. 12Hr
#6583 (Ca^{2+})	Tap water	Approx. 3Hr

6. Measuring

● Preliminary checks

- ① Remove the protective cap.
- ② If the internal solution level is low, remove the rubber stopper for the internal solution filling opening and add more internal solution to the outer tube.

● The measurement method

Create a calibration curve for the electrode using 2 or 3 types of standard solutions with 10-fold differences in concentration. This step is both for preparing the calibration curve, and also for verifying operation. The next step is to make quantitative measurements of the sample solutions using one of the following methods.

- (a) Direct Potentiometric Method
- (b) Known Addition Method
- (c) Known Subtraction Method
- (d) Potentiometric Titration Method

Consult the JIS K 0122 standards (General Rules for Ion Electrode Methods) which covers each of these measurement techniques in detail. Here is a brief overview of (a) the Direct Potentiometric Method.

With the Direct Potentiometric Method, you first create a calibration curve using two types of standard solutions with ion concentrations slightly higher and slightly lower than the ion concentration in the sample solution. You then measure the electric potential in the sample solution, and arrive at the quantification value using the calibration curve.

● Measurement operations (Also refer to the instruction manual for the tip type ion electrode)

- ① Wash the electrode well using pure (deionized) water, and wipe away any water droplets using gauze or tissue paper.
- ② Stir the solution at a constant speed (approximately 300 to 500 rpm) using the magnetic stirrer.

Note: If there is no supporting electrolyte in the test sample, or the concentration of the supporting electrolyte is too low, stirring the mixture too strongly will cause electric potential differences between the liquids and may make the electrode potential of the reference electrode unstable.

- ③ Immerse the electrode in the sample liquid. Be sure to keep the surface of the internal solution of the electrode above the surface level of the sample liquid.

Note:

- When taking measurements, remove the rubber stopper from the internal solution filling opening of the electrode and leave the internal solution cavity open to the atmosphere.
- The detection characteristic of the electrode varies with the temperature. You must maintain the standard liquid and the sample liquid at the same temperature during measurement.

● After measurements are completed

- ① Wash the electrode well with pure (deionized) water to remove any sample liquid, and completely remove any water droplets with gauze or tissue paper.
- ② Firmly plug the internal solution filling opening with the rubber stopper.
- ③ Attach the protective cap. (Make sure that the protective cap is completely dry.)

7. Storage

● Short-term storage

If you will be using the electrode again the next day, store the electrode by soaking it in the conditioning liquid used in step 5.(9). Store solid state membrane electrodes #8560 (Cl⁻) and #6561 (F⁻) by soaking them in pure (deionized) water.

● Long-term storage

- ① Remove the tip electrode from the combination type electrode and attach the rubber cap.
- ② Attach the protective cap. (Make sure the protective cap is completely dry.)
- ③ Store the tip type electrode together with the combination type electrode in a dry state.
- ④ The next time you use the electrode, start from the procedures given in section 5. "Before Usage".

8. Regular Maintenance

- ① If the electrode is used for long periods, the internal solution in the outer tube of the reference electrode may be contaminated by the sample liquid, or become diluted. Replace the internal solution for the reference electrode using steps 5.(2) to (3) once every one to five weeks.
- ② The ideal outflow of the internal solution is a slight oozing from the liquid junction. If the outflow is excessively small, then the electrode potential of the reference electrode will become unstable and will be more greatly affected by the stirring. When this happens, use steps 5.(4) to (7) to increase the outflow of the internal solution from the liquid junction. (Figure 5.)

Note: Make sure that the rubber cap is attached for the tip type ion electrode section. If this operation is performed with the tip electrode attached, it may cause damage to the electrode.

9. Supplementing the Internal Solution (Inner Tube)

Because the inner tube is sealed there is very little reduction in the amount of internal solution in the inner tube. For short-term usage of several months, there is no need to supplement the internal solution. However, for long-term usage, the internal solution will need to be supplemented when the amount of liquid in the inner tube is reduced to approximately half. (The frequency of replacement differs depending on the usage and storage conditions, however, when used at normal temperatures, replacement is required approximately once a year. When supplementing the internal solution in the inner tube, it is necessary to disassemble the electrode which will require special care.)

● Supplement method

(Disassembly)

- ① Remove the protective pipe and tip electrode from the combination type electrode, and attach the rubber cap to the tip electrode part.
- ② Remove the rubber stopper in the internal solution (outer tube) filling opening, and use the syringe to extract the internal solution.
- ③ Twist the electrode cap with your hand and move it 5 to 10 cm toward the electrode connector. (Figure 3).
- ④ Holding the electrode cap with your hand, stand the electrode vertically on a desk or other horizontal surface and press it down two or three times to release Body 1 from the inside of the electrode.
- ⑤ Move the silicon tube of Body 1 down to expose the internal solution filling opening (inner tube).
- ⑥ Using the syringe, fill the tube with gel type internal solution (#330) almost to the internal solution filling opening (inner tube).

(Assembly)

- ① Return the silicone tube to its original position and close the internal solution filling opening (inner tube). (Check that the opening is completely sealed.)
- ② If body 1, body 2 and the liquid junction are dirty, wash and clean them with pure (deionized) water.
- ③ Insert Body 1 into Body 2. (Check that the O-ring is firmly sealed.)
- ④ Return the spring to the top part of Body 1.
- ⑤ Holding the rubber sleeve with your hand, twist the electrode cap approximately 90 degrees, and remove the rubber sleeve from the electrode cap section. (Figure 4).
- ⑥ Arrange the electrode cap so that the HORIBA logo is facing in the same direction as the internal solution filling opening (outer tube), and attach the cap to Body 2.
- ⑦ Hold the electrode cable (connector side) in your hand, and move the rubber sleeve close to the electrode cap. After positioning the rubber sleeve so that it will fit in the hole of the electrode cap, attach it to the electrode cap.
- ⑧ Fix the rubber sleeve in place by twisting it approximately 90 degrees.
- ⑨ Using the syringe, supplement the internal solution manufactured in step 4. (1).
- ⑩ Follow step 5. (6) to start the outflow of the internal solution from the liquid junction.
- ⑪ Store according to the steps given in item 7. "Short-term storage".


Figure 2. Assembly and disassembly diagram


Figure 3. Assembly step (3)


Figure 4. Assembly step (5)


Figure 5. Starting the outflow of the internal solution