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
## ISFET pH sensor 0131 Instruction Manual

Thank you for purchasing a HORIBA ISFET pH sensor. ISFET sensor 0131 is applicable as a consumable part of ISFET Electrode 0030-10D. Read operation manual of ISFET pH sensor and this manual carefully before use.


### Cautions on use

- Handle carefully the sharp electrode tip to not get hurt.
- Do not scratch or apply excessive force to the sensing part and liquid junction.
- Do not touch the sensing part and sensor terminal. This may cause malfunction.
- Do not allow the electrode connector and sensor terminal to come in contact with water or wet hands.
- Before replacing the sensor or batteries, turn off the power and wipe water off the electrode.
- For the first use, or use after a long interval (one week or longer), immerse the sensor tip in 3.33 mol/L KCl solution (#300) for 30 minutes or longer before use.
- Make sure to use the electrode at a temperature within the specified usable range (sensor: 0 °C to 60 °C, converter: 0 °C to 40 °C).
- The product should not be subject to rapid temperature change.
- Dirt on the sensing part or liquid junction will affect accuracy adversely. Store the product in dry conditions. Storage area should be kept away from high temperature, high humidity, and direct sunlight.
- Do not use organic solvents such as acetone, tetrahydrofuran (THF) for cleaning. The sensor materials include resins, including ABS resin, epoxy resin, and polyethylene resin.
- White powder or solution may be generated at the liquid junction. This is leaked internal solution of the reference electrode, and not a symptom of trouble. Rinse the sensor with ion-exchange water before use.
- The sensor is a consumable. If it is damaged or deteriorated, replace a new sensor. (The sensor cannot be repaired).

### Packaged Contents

- Sensor (0131): 1
- Instruction manual (this book): 1

### Components

### Specifications


Sensor model	0131
Measurement range	pH 0 to 14
Usable temperature range	[Sensor] 0 °C to 60 °C
Storage temperature range	0 °C to 40 °C
Length	100 mm
Sensing device	ISFET (Ion-sensitive field-effect transistor)
Internal electrode	Silver/ silver chloride electrode
Reference electrode internal solution	Carboxymethyl cellulose (CMC), glycerin, supersaturated KCl
Liquid junction material	porous sintered polyethylene
Sensor material	Tantalum pentoxide, platinum, acrylonitrile-butadiene-styrene (ABS) resin, epoxy resin, polyethylene resin
Height of sensing part:	16 mm (from the tip of electrode)
Application	0030-10D

### Assembling sensor and converter


#### Note

**Do not touch the sensor terminal. It may be damaged by static electricity.**


#### 1 Turn OFF the power.


#### 2 Turn the nut to the left to remove it. And pull out the old sensor from the converter.

#### 3 Remove the end cap from the sensor, and the seal from the liquid junction. Make sure that the O-ring is attached to the sensor.

#### 4 Insert the sensor to the converter slowly checking that the O-ring is untwisted, and match the Δ marks on the sensor and converter. Then turn the nut to the right to fix the sensor.

#### 5 Immerse the sensing part and liquid junction completely into 3.33 mol/L KCl solution (#300) for 30 minutes or longer.