

Low Temperature Macro Photoluminescence and PLE System



- Modularity and versatility
- Many configurations available
- Fully automated devices
- Multichannel detection
- Controlled by *DataMax*[™]

Figure 1: *FluoroLog-3*[®] system incorporating a closed cycle helium cryostat

The *Jobin Yvon FluoroLog-3*[®] system makes an ideal turnkey macro PL / PLE system. The large format sample chamber allows the easy installation of a standard helium cryostat for low temperature studies.


The system is controlled by the *Jobin Yvon DataMax*[™] software, which can perform excitation, emission, time base and synchronous scans in both eV and nm units.

Instrument's Configuration

Excitation is from a powerful 450 W xenon lamp coupled to a double grating fast scanning monochromator. The double grating arrangement greatly reduces stray light and as the gratings are kinematically mounted they can easily be changed by the user to maximize efficiency at the desired wavelength.


The sample chamber optics, like the rest of the *FluoroLog-3*[®], consist only of mirrors which unlike lenses will maintain their efficiency and focus throughout a wide wavelength range.

The emission goes through a *Triax 552* 0.5 m spectrograph, which offers high resolution and can accommodate a variety of different detectors such as, photon counting PMT, CCD, ICCD, Ge diode and InGaAs arrays. The *Triax 552* is unique in that two array detectors (CCD and InGaAs) can be accommodated simultaneously.



Results

PL and PLE spectra of a GaN sample at various temperatures.



Sample courtesy of University of Limerick

Specifications

Source – 450 W xenon lamp

Excitation range – 230 nm to 1000 nm

Emission range - 230 nm to 3000 nm

“Need more information on these devices?
Contact us!”

In France :
Jobin Yvon S.A.
16-18, rue du Canal
91165 Longjumeau cedex
Tel: +33 (0)1 64 54 13 00
Fax: +33 (0)1 69 09 93 19
Email: fluorescence@jyhoriba.fr

In the USA :
Jobin Yvon Inc.
3880 Park Avenue
Edison, NJ 08820
Tel: +1 732 494 8660
Fax: +1 732 549 5157
Email: fluorescence@jyhoriba.com

China: +86 10 6836 6542
Germany: +49 89 46 23 17 0
Italy : +39 2 57 60 30 50
Japan: +81 35667 7351
U.K. : +44 20 8204 81 42

