

Only with Lite^{DM} - Now with multi-system connectivity



Connections* available include:

HORIBA Analyzers

- Micros 60
- Micros ES60
- Pentra 60 C+
- Pentra XL80
- Pentra C400

3rd Party Analyzers

- Siemens Clinitek Status[®]
- Abaxis Piccolo Xpress[™]
- Tosoh AIA 360[®]
- Quidel Triage[®]
- Various urine strip readers

Test and Test Panel Customization

- Define multiple panels
- Create qualitative tests
- Enter any off-line test results manually
- Reference ranges for gender and age

Simple Navigation and Reporting

- Simple and intuitive to learn and use
- CLIA compliant patient reports
- Single and cumulative patient reporting
- Quick and easy patient result history view
- Consolidated patient report with all results
- Add off-line point-of-care results
Flu | Strep | HCG | and others as needed
- Preliminary reports
- Easy backup and restore
- No monthly or annual maintenance fees

Additional Connectivity

- Connect directly to LIS systems
- Receive orders and transmit results to EMR systems

Connect up to 4 Analyzers

(including Micros 60) with one Lite^{DM}



Micros 60/ES 60



Pentra C400



*Possible interface connections

RESULT REPORT

HORIBA Medical
 9755 Research Drive
 Irvine CA 92618

DIRECTOR Edward Forrest, M.D.
CLIA 123D45678
Phone 306-987-5695
Fax 306-877-9023

Patient Palmer, Mary

Medical #: MI2345

Sample ID 111201
Doctor Patel, John

Age 28 **DOB** 10/12/1990


01/01/1999

Verified By 08/15/2016 5:15:28 PM
 USER


Gender F **Child**

	Test	Result	Abnormal	Flags	Normals	Units	Notes
Micros 60 Hematology Results	WBC		12.2	H	(3.5 - 10)	10 ³ /mm ³	
	RBC	5.19			(3.8 - 5.8)	10 ⁶ /mm ³	
	HGB	15.3			(11 - 16.5)	g/dl	
	HCT	45.6			(35 - 50)	%	
	MCV	88			(80 - 97)	µm ³	
	MCH	33.2			(26.5 - 33.5)	pg	
	MCHC	33.4			(31.5 - 35)	g/dl	
	RDW	12.4			(10 - 15)	%	
	PLT	281			(150 - 450)	10 ³ /mm ³	
	MPV	9.0			(6.5 - 11)	µm ³	
	LYM%	21.7			(17 - 48)	%	
	LYM#	2.60			(1.2 - 3.2)	10 ³ /mm ³	
	MON%		3.5	L	(4 - 6)	%	
	MON#	0.40			(0.3 - 0.8)	10 ³ /mm ³	
	GRA%	74.8			(43 - 76)	%	
GRA#		9.10	H	(1.2 - 6.8)	10 ³ /mm ³		
Accustrip Urinalysis Results	LEU		75	*	(-)	Leu/µl	
	NIT		POS	*	(-)	mg/dl	
	PRO		100	*	(-)	mg/dl	
	pH	6, 5			(-)		
	BLO	NEG			(-)	Ery/µl	
	SG	1.015			(-)		
	KET		25	*	(-)	mg/dl	
	BIL		4	*	(-)	mg/dl	
	GLU		50	*	(-)	mg/dl	
	UBG		8	*	(-)	mg/dl	
Manual Entry Test	PREG	NEG			(-)		ww


WBC Histogram



RBC Histogram



PLI Histogram



Flags
 WBC Flags L1 M2G1 G3
 PLT Flags SCL
Notes

Micros 60 Specimen Type : whole blood
 Micros 60 Analyzed D/T : 07/27/2008 05:14 PM
 Micros 60 Collected D/T : 08/15/2016 05:01 PM
 Micros 60 Received D/T : 08/15/2016 05:01 PM

Reviewed By: _____ Date: _____

MED-ADV-10244 V1.0