


## SPECTRAL EFFICIENCY CURVE

### Reference 510 91 ---

Grating Type	Plane Mechanically Ruled Grating
Groove Density (gr/mm)	900
Spectral Range (nm)	
Blaze Wavelength (nm)	350
Blaze Angle (deg)	9°00'
Configuration	Quasi-littrow

### 510 91 ---



This efficiency curve should be used as reference only and does not imply any guarantee on actual performance. Actual performance of the grating will depend on several variables including the usage geometry and the coating. Efficiency curves labeled as "Theoretical" are calculated estimates of performance based on groove profile.